

## Western Immunoblotting Protocol

### **Blocking**

1. Block membrane by incubating 1 hour at room temperature or overnight at 4°C with shaking in Blocking Solution (5% BSA or nonfat milk in TBST (50mM Tris, 100mM NaCl, 0.05% Tween-20, pH 7.6)). Note: Use 5% BSA in Blocking Solution for phosphor specific antibodies.

### **Incubation with Primary Antibody**

2. Dilute primary antibody at the appropriate dilution in Blocking Solution.  
3. Incubate the membrane with diluted primary antibody for 1 hour at 37°C, or 2 hours at room temperature, or overnight at 4°C with agitation.  
4. Remove antibody solution. Wash the membrane 3 times for 5-10 minutes each time at room temperature in TBST (50mM Tris, 100mM NaCl, 0.05% Tween-20, pH 7.6) with shaking. Note: Increase the concentration of Tween-20 to 0.1% reduces the background and increases the specificity, but it will reduce the sensitivity.

### **Incubation with Second Antibody**

5. Incubate membrane with secondary AP or HRP conjugate diluted (according to manufacturer's instructions) in Blocking Solution for 1 hour at room temperature with shaking.  
6. Wash the membrane as Step 4.  
7. Wash membrane with TBS for 2-5 minutes before proceeding Chemiluminescent Reaction.

### **Chemiluminescent Reaction**

8. Prepare and use the Chemiluminescent substrate (for AP or HRP) according to the manufacturer's instructions.  
9. Immediately wrap the membrane and expose to X-ray films for 10 second to 1 hour period. The exposure time may vary according to the amount of antibody and antigen.

### **Peptide Competition**

Before proceeding Western Immunoblotting, add Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.